

FAQ: “Can HunterLab EasyMatch QC and Universal software programs be installed on the same hard drive? And of course, be operated by different people?”

Yes, both our current EasyMatch QC and legacy Universal software packages can be run on the same PC but not at the same time as Universal does not like to share communications with the same sensor. Just close one program before starting the other, and there will be no conflicts.

The PC should have a WIN XP operating system to operate both color quality packages simultaneously.

Universe was qualified for operation only up to WIN XP when it was replaced by EasyMatch QC.

EZMQC will operate on WIN XP and all Windows [current operating systems](#).