

EZMQC Program Launch RunTime Error

I keep getting a Runtime Error message whenever I try to launch the program.

The Error Tool Box Says "Microsoft Visual C++ Runtime Library." Runtime Error!

C:/Program Files > HunterLab > EasyMatchQC > EasyMatchQC.exe

This application has requested the Runtime Error to terminate it in an unusual way. Please contact the application's support team for more information.

Use the EZMQC Installation CD and when the installation screen appears select "Repair" then select next. This should reload the run-time libraries needed to operate the software.